



DURAN® Rundkolben, enghalsig, mit Bördelrand

DURAN® Round bottom flasks, narrow neck, with beaded rim

	Inhalt [ml] Capacity [ml]	d [mm]	d ₁ [mm]	h [mm]	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
DWK	50	51	26	95	10	21721170
	100	64	26	110	10	21721240
	250	85	34	144	10	21721360
	500	105	34	168	10	21721440
	1000	131	42	200	10	21721540
	2000	166	42	260	10	21721640 ¹
	3000	185	50	260	1	21721680 ¹
	4000	207	52	290	1	21721710
	5000	223	50	305	1	21721730 ¹
	6000	236	51	355	1	21721770 ¹
	10000	279	65	380	1	21721860
	12000	295	65	380	1	21721870 ^{1,2}
	20000	345	76	515	1	21721910 ²
	50000	480	110	600	1	21792157 ³

- Nach DIN ISO 1773
- Rückverfolgungscode
- USP Standard
- ¹ Nicht nach DIN ISO
- ² Nach ASTM E 1403
- ³ ohne Druck
- Temperaturbeständigkeit: + 200°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: Borosilikatglas 3.3

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- According to DIN ISO 1773
- Retrace code
- USP Standard
- ¹ Non-DIN ISO size
- ² Conforms to ASTM E 1403
- ³ without print
- Temperature resistance: + 200°C
- Pressure resistance: -1/+0.5 barg
- Material in contact with product: Borosilicate glass 3.3

Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.