



DURAN® Rundkolben, weithalsig, mit Bördelrand

DURAN® Round bottom flasks, wide neck, with beaded rim

	Inhalt [ml] Capacity [ml]	d [mm]	d ₁ [mm]	h [mm]	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
DWK	50	51	34	105	10	21741170 ¹
	100	64	35	110	10	21741240
	250	85	51	143	10	21741360
	500	105	50	168	10	21741440
	1000	131	50	200	10	21741540
	1000	131	65	200	10	21741550 ¹
	2000	165	76	240	10	21741630
	2000	166	50	240	10	21741640 ¹
	3000	185	65	260	1	21741680 ¹
	4000	206	76	290	1	21741710
	5000	223	65	310	1	21741730 ¹
	6000	236	89	330	1	21741760
	6000	236	65	330	1	21741770 ¹
	10000	279	89	420	1	21741860 ¹
	20000	345	89	520	1	21741910 ¹

- Nach DIN ISO 24450
- USP Standard
- ¹ Nicht nach DIN ISO.
- Temperaturbeständigkeit: + 200°C
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: Borosilikatglas 3.3

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- According to DIN ISO 24450
- USP Standard
- ¹ Non-DIN ISO size
- Temperature resistance: + 200°C
- Pressure resistance: -1/+0.5 barg
- Material in contact with product: Borosilicate glass 3.3

Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.