



NS-DURAN® Gewinderohr mit GL-Gewinde und NS-Kern

NS-DURAN® GL-Screwthread to ground joint adapter

	DIN-Gewinde GL DIN thread GL	L [mm]	Kern NS Cone size NS	Verpackungseinheit/Stück Pack/Quantity	Artikelnr. Article no.
DWK	14	30	14/23	10	24840620
	14	35	19/26	10	24840720
	14	40	24/29	10	24840820 ¹
	14	40	29/32	10	24840920 ¹
	18	35	14/23	10	24841610 ¹
	18	35	19/26	10	24841710
	18	40	24/29	10	24841810 ¹
	18	40	29/32	10	24841910
	25	40	19/26	10	24842720 ¹
	25	40	24/29	10	24842820 ¹
	25	40	29/32	10	24842920
	32	50	24/29	10	24844810 ¹
	32	50	29/32	10	24844910

- DWK geschliffen - 600er-Körnung
- Schliff fett benötigt - alternativ PTFE-Dichtung, S. 43/44
- Temperaturbeständigkeit: + 200°C - nicht isoliert, Minustemperatur limitiert durch Schliff fett/dichtung
- Druckbeständigkeit: -1/+0,1 barg bei Verwendung herkömmlicher Klemmen, -1/+0,5 barg mit NS-Schliffverschraubung, S. 38
- Material produktberührend: Borosilikatglas 3.3
- ¹ Nicht nach DIN

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- DWK ground - 600 grain
- Ground joint grease required - alternatively PTFE seal, p. 43/44
- Temperature resistance: + 200°C - not insulated, minus temperature limited by grinding grease/sealing
- Pressure resistance: -1/+0.1 barg when using conventional clamps, -1/+0.5 barg with NS ground joint p. 38
- Material in contact with product: Borosilicate glass 3.3
- ¹ Non-DIN size

Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.