



NORMAG



VH-Durchgangs-Spindel- ventile mit Temperiermantel

VH straight spindle valves with tempering jacket

	Nenngröße NS Nominal size NS	Rohr d1 x s1 [mm]	Rohr D x S [mm]	Artikelnr. komplett Article no. complete	Artikelnr. Glasteil Article no. glass part
NORMAG	3	24 x 1,8	8 x 1,5	VHTS032418	VHTS032418G
	6	28 x 2,0	10 x 1,5	VHTS062820	VHTS062820G
	10	34 x 2,0	15 x 1,8	VHTS103420	VHTS103420G
	15	48 x 2,3	22 x 2,5	VHTS154823	VHTS154823G

VH-Ventilstempel für Ventile mit Temperiermantel VH valve stamps for valves with tempering jacket					
	Nenngröße NS Nominal size NS	d1 [mm]	d2 [mm]	L [mm]	Artikelnr. Article no.
NORMAG	3	8	7	36	VHT03S
	6	10	9	49	VHT06S
	10	16	15	62	VHT10S
	15	22	21	75	VHT15S

- Weltweit niedrigste Leckrate von 1×10^{-12} mbar*l/s
- Totraumfrei optimiert
- Flexibel einsetzbar für Flüssigkeiten und Gase
- Sicherheit durch Ausdrehschutz
- Durch komplexes aber unauffälliges Dichtungssystem unter PTFE-Überzug hochvakuumdicht, ohne offene Dichtringe, fettfrei und totraumfrei optimiert
- Dichtsystem nachstellbar
- Nur noch in starkwandig und extra starkwandig erhältlich
- Material Handgriff nicht produktberührend
- Material Handgriff: PA 6 GF 30, 0/+100°C
- Material O-Ring: FKM, -40/+180°C
- Material Metallteile: Messing, -50/+200°C
- Temperaturbeständigkeit: -100/+200°C, bei T < -50°C kann ein nachstellen notwendig sein
- Druckbeständigkeit: -1/+0,5 barg
- Material produktberührend: Borosilikatglas 3.3, PTFE
- Mit Temperiermantel

Hinweis: Für Glasartikel, die durch Weiterverarbeitung von der Originalform abweichen, sind den technischen Regelwerken und Vorschriften entsprechend, die zulässigen Betriebsdaten zu ermitteln.

- Worldwide lowest leakage rate of 1×10^{-12} mbar*l/s
- Optimized without dead space
- Flexible use for liquids and gases
- Safety through protection against unscrewing
- Highly vacuum-tight, without open sealing rings, grease-free and free of dead space due to complex but inconspicuous sealing system under PTFE coating
- Sealing system adjustable
- Only available in thick-walled and extra thick-walled
- Material handle not in contact with product
- Material handle: PA 6 GF 30, 0/+100°C
- Material O-ring: FKM, -40/+180°C
- Material Metal parts: brass, -50/+200°C
- Temperature resistance: -100/+200°C, at T < -50°C readjustment may be necessary
- Pressure resistance: -1/+0,5 barg
- Material in contact with product: Borosilicate glass 3.3, PTFE
- With tempering jacket

Important: For glass articles which differ from the original shape after further processing, the permissible data must be determined in accordance with the technical rules and regulations.